**PUSH ORTHO THUMB BRACE CMC**

**INTRODUCTION**  
Push ortho products offer effective solutions for complaints or reduced functioning of the locomotor apparatus. The products have been approved by medical specialists and scientists.

Your new thumb brace is a Push ortho product, characterised by a high level of effectiveness, comfort and quality. The Push thumb brace CMC guarantees improved stability and pain relief. It increases the ability to bear physical stress. Its optimal fit provides maximum mechanical support.

Push ortho products are innovative and are carefully manufactured for you. This guarantees maximum security and greater freedom of movement.

Read the product information carefully before you use the Push CMC. If you have any questions, please consult your reseller or send an e-mail to [info@push.eu](mailto:info@push.eu).

**PRODUCT CHARACTERISTICS**  
Osteoarthritis of the carpometacarpal joint of the thumb (CMC-1) causes functional instability of the joint at the base of the thumb as a result of a weakening of the ligaments.

The Push CMC stabilises the basal joint of the thumb (CMC-1) and places the thumb's metacarpal (metacarpal 1) in a functional position.

This creates a functional thumb arc, which remains stable while the hand is used. It relieves the pain when the thumb is under stress. This method of stabilisation leaves the MCP-1 joint and the wrist free, to allow maximum functioning of the hand. For optimal performance, the aluminium core needs to fit closely to the thumb. The specialist supplying the product can, if necessary, adjust the brace individually by reshaping the aluminium core. The selected materials and rounded edges guarantee maximum comfort. The simple closing system enables you to fit the brace yourself without any effort.

The Push CMC is a slim and compact product. The brace can be used in the water. After use, simply let it dry before refitting it. As hygiene is important around the hand, the synthetic material is antibacterial and the brace can be machine washed at 40°C.

- INDICATIONS**
- CMC-1 osteoarthritis
  - Postoperatively after arthroplasty of the CMC-1 joint
  - Instability of the CMC-1 joint

There are no known disorders where the brace cannot or should not be worn. In the case of poor blood circulation in the limbs, a brace may cause hyperaemia. Avoid applying the brace too tightly.

Pain is always a sign of warning. If the pain persists or the complaint deteriorates, it is advisable to consult your doctor or therapist.

**FITTING**  
The Velcro fasteners are marked with symbols. These symbols will help you close the Velcro fastenings in the proper order. First Symbol  , then Symbol .

- 1 Slide the brace around your thumb, with the bands at the back of the hand. The aluminium core then encloses the base of the thumb. The Velcro of the bands must be loose, but the bands may remain in the holes.
- 2 Keep your fingers and thumb bent in a functional position (see Figure B) and close the lower band nearest to the wrist. (Symbol  )
- 3 Close the upper band (Symbol  ).

**CARE AND MAINTENANCE**  

Before washing, close the Velcro to protect other laundry.

The brace can be washed at temperatures up to 40°C, on a delicate fabrics cycle or by hand. Do not use any bleaching agents.

A gentle spin-dry cycle can be used: afterwards, dry the brace naturally. (Do not dry it on a radiator or in a tumble dryer).

To extend the brace's service life, it should be dried after use, if necessary, before being stored.

**IMPORTANT**  
Correct use of your brace requires a detailed diagnosis of the cause of your pain complaints or instability. Always consult your GP, a physiotherapist or a specialist before applying the brace.

- For the brace to function properly, the following requirements must be met:
- The brace fits properly. If necessary, an acknowledged expert should adapt the brace to fit.
  - The washing instructions must be followed
  - The Velcro fasteners must still work
  - There is no damage (tears, disconnected components)

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