## PUSH ORTHO ANKLE FOOT ORTHOSIS AFO

### INTRODUCTION

Push ortho products provide effective solutions in the case of a disorder or reduced functioning of the locomotor apparatus. Approved by physicians and scientists.

Your new ankle-foot orthosis is a Push ortho product, characterised by a high level of effectiveness, comfort and quality.

Push ortho is used for complaints or reduced function of the locomotor apparatus. Adequate mechanical support is key. Push ortho can play an important role in cases of repositioning and stabilization thus compensating for loss of normal function.

Push ortho products are innovative and are manufactured for you with great care. This ensures maximum security and more freedom of movement. To promote satisfied use.

Read the product information carefully before you use the Push AFO. If you have any questions, please consult your reseller or send an e-mail to [info@push.eu](mailto:info@push.eu).

### PRODUCT CHARACTERISTICS

The Push AFO provides effective support in the case of foot drop, ensuring minimal disturbance of the person's natural gait, thus allowing safe and efficient walking.

The Push AFO is an ankle foot orthosis that lifts the foot in a unique combination of functions during the swing phase, slows down the plantar flexion after heel contact and subsequently releases the dorsal flexion of the ankle joint during completion of the step. The orthosis also provides mediolateral support for the ankle.

The functional strap system is continuously adjustable, allowing the degree of dynamic foot lifting during the swing phase to be adapted to the user's individual needs. During the landing phase, the plantar flexion is controlled by the braking function of the elastic straps. The relatively flexible sole of the Push AFO allows dorsiflexion in the ankle joint and makes it possible to use the residual function of the calf muscles during the completion phase of the step. This creates an efficient and natural gait.

The Push ortho works in different types of shoes and does not require a special (last-made) shoe. Whether the shoe usually provides fixation of the orthosis around the foot, the Push AFO itself possesses a fitting strap system. This innovative construction allows the Push AFO to be used even with slippers.

The fibreglass-reinforced AFO is a lightweight construction that provides freedom for the resting function of the muscles of the lower leg, as a result of which the product is experienced as very comfortable.

At night the Push AFO can be used as a night splint to prevent cramping or hanging of the foot. It is advised not to overtighten the non-elastic strap system and not to close the elastic straps crosswise but to close them loosely along the side of the leg.

### INDICATIONS

- Loss of control of the ankle dorsiflexors during the swing phase and during loading response (heelstrike / foot flat) as a result of failure of the tibialis anterior, possibly in conjunction with weakness of the tibialis posterior.
- Disorders possibly involving loss of the function of the ankle dorsiflexors are: CVA, hernia, multiple sclerosis, postpolio syndrome, muscular dystrophy, spinal cord injury, polyuropathy.
- **Note:** in the case of use by diabetics, a thorough assessment by a care professional of potential pressure points is required.

### CONTRAINDICATIONS

- Loss of lack of control of the calf muscles, moderate to severe spasticity, extreme pes equinus, foot deformity.

### FITTING

- Open all straps (A).
- Insert your foot underneath the crossed straps (A).
- Place your heel securely at the back of the heel area (B).
- Close the upper closure (around the lower leg) (C).
- Apply the correct tension to the inelastic straps (C), pull the closure (C) upwards and fix it at the back of the leg.
- Close the elastic strap from the outside of the foot to the inside of the calf (D), ensuring a comfortable tension (C).
- Then close the second elastic strap (C).
- You may need to retighten and close the inelastic straps (C) at this stage.

### INSTRUCTIONS FOR INDIVIDUAL FITTING BY A PROFESSIONAL SPECIALIST

#### Determining size

The Push AFO is available in three sizes in a left and right execution. Place the foot at a 90-degree angle in relation to the lower leg, in order to measure the exact size. Then measure the circumference around the heel and the instep of the foot as shown in the image. See table for the correct size at each foot circumference.

To be on the safe side, also check whether the measured size matches the user's shoe size.

#### Customizing

The Push AFO must at various points be adjusted to size by a fitting professional.

#### The sole

Adjust the sole of the frame in such a way that it ends between the ball of the foot and the toes, without extending beyond the sides of the foot. The softshell sole is also customizable. Keep it just slightly longer than the frame sole, so that the edge of the frame is well padded.


#### The front closure

To adjust the front closure to the right size, the removable end of the front closure (1) should be taken off the front closure. Next the band must be cut to size (2) with a pair of scissors and the end must be put back again (3).

#### Individual setting

- To prevent users with sensory loss over-tightening the straps, the attached setting indicators are to be applied to the loop strap at the back and at both sides of the orthosis (4). In this way the individual setting can easily be found every day.
- In particular in the case of use by diabetics, the creation of folds or other pressure points must be avoided during fitting.

#### WASHING INSTRUCTIONS

Before washing, close the Velcro to protect the brace and other laundry. The brace can be washed at temperatures up to 30 °C, on a delicate fabrics cycle or by hand. Do not use any bleaching agents.

A gentle spin-dry cycle can be used: after this, hang the brace up to dry in the open air (not near a heating source or in a tumble drier).

To extend the brace's service life, it should be dried after use, if necessary, before being stored.

#### IMPORTANT

The efficient use of Push AFO depends on the correct diagnosis. Always check first with a specialist before using a orthosis.

For the orthosis to function properly, please note the following:

- special attention needs to be given to individuals suffering from loss of sensitivity
- the Push AFO is worn inside a shoe
- to follow the washing instructions
- that the elastic material can still be stretched
- the velcro fastens efficiently
- the orthosis is not severely damaged (no tears, loose seams)
- the orthosis is not subjected to temperatures higher than 70 °C
- the use of a rigid orthosis within a shoe may cause accelerated wear of the inside of the shoe.

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